



## HellabioScan Software v.03 Film Analysis Quick Guide

# Hemoglobin electrophoresis film analysis

### Scan Films

HellabioScan analyzes images that are stored at any picture format (\*.jpg, \*.pif, \*.tif, \*.bmp etc) in your computer. The user can use any peripheral device in order to obtain these images (such as Scanner, Digital Cameras etc.) and store them to his computer.

Recommended analysis for a typical Electrophoresis film is **100 dpi** by 256 gray level color depth (it combines a very good final result with a small image size in computer's hard disk). At this point we must refer that the basic scanning parameters such as dpi analysis or color depth have small effect to the final results of HellabioScan software.

### Saving Images

After scanning completion HellabioScan program asks the user to save the image. Recommended storage file → Follow the path **C: →program files (X86)→evresis→HellabioScan(SQL) → Images**. Alternatively the user can save the archives to a folder of his choice.

After saving completion HellabioScan captures automatically the scanned image and presents it on the Application's Desktop (default position is on the right side of desktop). **Figure 1** shows HellabioScan desktop after a typical electrophoresis film scanning.





Fig. 1

Additionally the user can load a saved image at any time, using the Open Image selection in File submenu (Ctrl + O) or by click on Toolbar's Open Image  Icon.


### **Importing patient data**

After scanning and saving completion, the user should enter the patient's personal data. HellabioScan software uses a database for patient registration. Using the option **tables** (above the toolbar), select the **Patients** option (**Fig. 2**)




**Fig. 2**

The database table will appear (**Fig. 3**).



**Fig. 3**

Using the option **Tables** again, select the **New Record** option, from the Tables submenu (**Fig. 4**)



**Fig. 4**

The program displays a submission form. The user should fill in all details, for each patient. **Figure 5** shows the submission form. Click on the option **Save**, to register the new patient.

A screenshot of a 'Patients' submission form. The form has a title bar with a person icon and the text 'Patients'. It contains several input fields: 'SSN', 'PersonName', 'Address', 'City', 'Country', 'Phone', 'Date of Birth' (with a dropdown menu showing '1/1/1990'), 'Doctor's Name', 'Laboratory', and 'Comments'. To the right of the 'Country' and 'Phone' fields is a button labeled 'Protocols' with a small icon. At the bottom of the form are two buttons: 'Save' (with a floppy disk icon) and 'Cancel' (with a red X icon).

**Fig. 5**

## Selecting Lanes

In general occasion Electrophoresis films contains more than one lane. The user can select for analyzing one or more lanes by drawing a rectangle using the computer's mouse. Just click on upper left corner and drag the mouse to the bottom right corner, drawing a rectangle that contains the lanes you want to analyze (**Fig. 6, 7, 8**).




Fig. 6




Fig. 7





Fig. 8

After the lane selection the user can analyze the selected lanes using the Analyze button  on Toolbar or choosing **Analyze** in File submenu (**Fig. 9**).

The program asks about the type of Electrophoresis film. Select **Hemoglobin (in alkali or in acid)** Electrophoresis Type (**Fig. 10**).




Fig. 9




Fig. 10

After the type selection the program automatically appears as many Calculation Windows as the selected lanes of the film. Figure 11 shows the HellabioScan Desktop after a 4-lanes selection of **4 blood samples**.





Fig. 11

Every Calculation window is an independent object. User can minimize, move or close any of the windows. Every Calculation window is labeled by a number according to the number of the lane in the film. In a Calculation Window the user can input several data about the examined sample. The **Patient** field must contain information about the examined sample (e.g. *the patient's name*) because it will be the filename of the specific examination. This is the name that the user must search in the future, if he wants to retrieve data about this examination.

Using the button  the user can enter the patient's name (Fig. 12).






Fig. 12


Select the patient by double-clicking on the arrow that appears on the left side of the table (**Fig. 13**).



ID	SSN	PersonName	Address	City	Country	Phone	BirthDay
45a423b9-		Sample 1					16/2/1997
7e10e7af-9		Sample 3					1/1/1990
898de7d1-		Sample 2					1/1/1990

**Fig. 13**


Using the Calculate button  on Toolbar the program executes an image processing algorithm that automatically defines the fractions of the image and calculates the density of each area. The results are presented at the down white part of the calculation window (**Fig. 14**).



Fractions	Values(%)	g/dl	NV(%)
HbA	97,60	0,00	96,7 - 98,5
HbA2	2,40	0,00	1,5 - 3,5

**Fig. 14**

The user can enter various data information's (e.g. hemoglobin value, comments concerning the diagnosis) (**Fig. 15**).



Patient: sample 1      Date: 12/10/2019

SSN: 012345

Type: Hemoglobin (in Alkali)      Total Quantity: 13.8

Comments: [Empty text box]


Sample Result Image: [Electropherogram image]

Buttons: PrintPreview, Print, Save

Fractions	Values(%)	g/dl	NV(%)
HbA	97,60	12,69	96,7 - 98,5
HbA2	2,40	0,31	1,5 - 3,5

**Fig. 15**

After the electropherogram analysis, the user can **Save** and **Print** the protocol (**Fig. 16**).



Patient: sample 1      Date: [Empty]

SSN: 012345

Type: Hemoglobin (in Alkali)      Total Quantity: 13.8

Comments: [Empty text box]

Sample Result Image: [Electropherogram image]

Buttons: PrintPreview, Print, Save (highlighted)

Fractions	Values(%)	g/dl	NV(%)
HbA	97,60	12,69	96,7 - 98,5
HbA2	2,40	0,31	1,5 - 3,5

**Fig. 16**




Fig. 16

### Adding an image to the protocol

Before printing, the user has the ability to preview the protocol (**Fig. 17**). HellabioScan software is automatically inserting the image of the electropherogram into the protocol. The user can choose whether or not to display the image in the protocol (**Fig. 18, 19**).




Fig. 17



Labels & Normal Values

Band Label	Normal Values	
HbA	96,7	98,5
HbA2	1,5	3,5

Header Image


Secondary Image

Fig. 18 The scanned film will appear in the protocol.

## Hemoglobin Electrophoresis (in Alkali)

Patient : sample 1  
SSN : 012345  
Date :

Examined Sample



Fractions	Values(%)	g/dL	NV(%)
HbA	97,60	13,47	96,7 - 98,5
HbA2	2,40	0,33	1,5 - 3,5

Comments:  
The Doctor

Copyright 2008-2015 EVRESIS S.A. Computer Software

Hemoglobin electrophoresis protocol

Labels & Normal Values

Band Label	Normal Values	
HbA	96,7	98,5
HbA2	1,5	3,5

Header Image


Secondary Image

Fig. 19 Blank window. The scanned image will not appear in the protocol.

## Hemoglobin Electrophoresis (in Alkali)

Patient : sample 1  
SSN : 012345  
Date :

Examined Sample



Fractions	Values(%)	g/dL	NV(%)
HbA	97,60	13,47	96,7 - 98,5
HbA2	2,40	0,33	1,5 - 3,5

Comments:  
The Doctor

Copyright 2008-2015 EVRESIS S.A. Computer Software

Hemoglobin electrophoresis protocol

## Hemoglobin variants

Hemoglobin variants are mutant forms of hemoglobin in a population (usually of humans), caused by variations in genetics. HellabioScan program performs an image processing algorithm that automatically defines the fractions of the image and calculates the density of each area. However, recognition of the hemoglobin variant is not supported by the software. User has to identify the type of hemoglobin variant and enter the data to create the protocol (**see example below**).

## Example

The screenshot shows the HellabioScan software interface. At the top, there are input fields for Patient (sample 1), SSN (012345), Date, Type (Hemoglobin (in Alkali)), and Total Quantity. Below these is a Comments field. In the center, there is a 'Sample' image showing a gel strip and a 'Result Image' showing an electropherogram with three peaks. To the right of the electropherogram are buttons for 'PrintPreview', 'Print', and 'Save'. Below the electropherogram is a table of results.

Fractions	Values(%)	g/dl	NV(%)
HbA	59,60	0,00	96,7 - 98,5
HbF/S	38,11	0,00	0,0 - 1,0
HbA2	2,29	0,00	1,5 - 3,5

Electropherogram of pathological sample. An extra band appears between HbA and HbA<sub>2</sub> (HbS, HbG, HbD, Lepore). After electrophoresis on acidic pH, the hemoglobin variant was identified as HbS.

## Protocol creation

After analyzing the electropherogram, the following data appear in the table of results (Fig. I)

Fractions	Values(%)	g/dl	NV(%)
HbA	59,60	0,00	96,7 - 98,5
HbF/S	38,11	0,00	0,0 - 1,0
HbA2	2,29	0,00	1,5 - 3,5

Fig. I

By clicking on the **PrintPreview** option, the user can modify the data (**Fig. II**). After modification, the user can **Save** and **Print** the protocol.

The main application window displays patient information and a chromatogram. The patient data includes:

- Patient: sample 1
- SSN: 012345
- Type: Hemoglobin (in Alkali)
- Total Quantity: (empty field)
- Comments: (empty text area)

The chromatogram shows a green shaded area under a curve with two peaks. The 'PrintPreview' button is highlighted with a red box.

The 'Labels & Normal Values' dialog box allows users to modify band labels and normal values. The 'HbF/S' band label is highlighted with a red box.

Band Label	Normal Values	
HbA	96.7	98.5
HbF/S	0.0	1.0
HbA2	1.5	3.5

Header Image: (empty field)  
Secondary Image: (empty field)

Buttons: Ok, Cancel

The 'Labels & Normal Values' dialog box shows the 'HbS' band label highlighted with a red box, and its normal values set to 0.0.

Band Label	Normal Values	
HbA	96.7	98.5
HbS	0.0	0.0
HbA2	1.5	3.5

Header Image: (empty field)  
Secondary Image: (empty field)

Buttons: Ok, Cancel

**Fig. II HbS Normal values NV (%): 0%**


# Hemoglobin Electrophoresis (in Alkali)

Patient : sample 1

SSN : 012345

Date :

Examined Sample



Fractions	Values(%)	g/dL	NV(%)
HbA	59,64	0,00	96,7 - 98,5
HbS	37,99	0,00	0,0 - 0,0
HbA2	2,37	0,00	1,5 - 3,5

Comments:

The Doctor



Copyright 2008-2015 EVRESIS S.A. Computer Software


## Hemoglobin electrophoresis protocol

### Manual correction

In some cases the software may not define the fractions properly. This may be due to technical errors (improper sample application, gel injury, insufficient staining-destaining). In this case the software allows the user to correct the errors manually.

### Toolbar

In the toolbar, the user has the following options:



Manual correction



Noise cleaning



Add a fraction separator




Undo last operation




Delete a fraction separator


### Example 1



**Displaced position of separator.**


For manually correction, the user must follow the following steps:

- I. Click on the option 
- II. Move the mouse to the displaced separator. A cross will appear on the separator (**Fig. 1.1**).




**Fig. 1.1**

- III. Move the separator (**Fig. 1.2**)



**Fig. 1.2**


## Example 2



### **More separator /separators**


For manually correction, the user must follow the following steps:

- I. Click on the option
- II. Click on the option
- III. The following table will display on the desktop (**Fig. 2.1**).




**Fig. 2.1**

- IV. Select the separator you want to delete (**Separator Number**). Click on the option **Delete**.
- V. Alternatively move the mouse over the separator you want to delete. A cross will appear on the separator.
- I. Right-clicking will display the following table (**Fig 2.2**)



**Fig. 2.2**

- II. Click on the option **Yes**. The separator has been deleted (**Fig. 2.3**).



**Fig. 2.3**

*If you need further assistance or have any questions related to the HellabioScan program, please do not hesitate to contact our team.*

**Email:** [hellabio@hellabio.com](mailto:hellabio@hellabio.com) &  
[hellabio@thermokoitida.gr](mailto:hellabio@thermokoitida.gr)