

# HEMOGLOBIN F (HbF) CONTROL

REF HbF01

LOT



HELLABIO Hemoglobin HbF01 Control is an assayed quality control material intended to be used in Hemoglobin electrophoresis in alkali and in acid. It is packaged in bottles, each containing 0.25 ml of control solution.

*Active Ingredients: HELLABIO HbF01 control contains a solution of unfixed stabilized human erythrocytes and preservatives.*

## PRECAUTIONS:

*This product contains human source materials that have tested non-reactive for hepatitis B surface antigen (HBsAg), hepatitis B core antibody (anti-HBc), hepatitis C virus (HCV), anti-HIV-1, anti-HIV-2, and human T-cell lymphotropic virus type 1 (HTLV-1).*

*Since no known test method can offer complete assurance that specimens of human origin will not transmit infectious disease, this control should be handled as a potentially biohazardous material. Follow the recommended "Universal Precautions" when handling this and any blood product.*

*Disposal of this product should be according to the practices of your institution for infectious medical waste. Discard all materials in a safe and acceptable manner and in compliance with all federal, state and local requirements.*

## STORAGE

The expiration date stated on the bottle is for product stored in freezing conditions (**under -8 °C**).

Bottles should be tightly closed after each use. Avoid exposing the bottle to temperatures greater than 8 °C. The solution is **ready to use**.

## INSTRUCTIONS FOR USE

For the proper use of Hemoglobin control, follow the steps below:

### For manual system:

- Place the vial on the laboratory workbench.
- Thoroughly mix by gently inverting the bottle and rolling between the palms until all cellular components are completely suspended. Do not shake or use a mechanical mixer.
- Twist open the cap.

**d)** Apply the control sample (**Ready to Use**) on the electrophoresis gel according to the instructions for use, (**user's Manual** kit HELLABIO HE, MHE, GHE, MGHE), for **30 seconds**.

**e)** Return the opened bottle to the chosen open-bottle storage temperature immediately after use. Do not use control solutions after the expiration date printed on the bottles and box.

### For HEELECS-1 instrument:

- Place the vial on the laboratory workbench.
- Thoroughly mix by gently inverting the bottle and rolling between the palms until all cellular components are completely suspended. Do not shake or use a mechanical mixer.
- Twist open the cap.
- Transfer 20µl of the sample control (**Ready to Use**) to one hole of the sample holder. Follow IFU (HEA13) for electrophoresis procedure.
- Return the opened bottle to the chosen open-bottle storage temperature immediately after use. Do not use control solutions after the expiration date printed on the bottles.

## LIMITATIONS

This product is intended for use as a quality control material in Hemoglobin electrophoresis in alkali and in acid. It is not for use as a calibration standard. Exposure to temperatures greater than 25 °C may affect product performance

## EXPECTED VALUES

100% HbF from umbilical cord.

## Hellabio

Prepared: Gkinis S.

Approved: Dr. Dimitriadis I.



DIMITRIADIS IOANNIS AND SON PC  
Production of electrophoresis products  
Business Incubator Thermi  
Steliou Kazantzidi 47, 57001, Thermi, Greece  
E-mail: hellabio@hellabio.com  
Phone: + 30 2311 999911 Fax: +30 2311 999912  
www.hellabio.com

Identification number: IFU/ HbF control

Updating no: 30 / 12-10-2022