



SAFETY DATA SHEET

Date of issue: 15.01.2025

According to Regulation (EC) No. 2020/878 & 1272/2008

Page: 1/7

SECTION 1: Identification of the substance/mixture and of the company

1.1 Product Identifier

REF: HCS01, HCS02

Product name: Human Serum Control

REACH registration number: The annual tonnage does not require registration

Kit Components/Mixtures (Buffers):

- Human serum control (1-vial)

1.2 Relevant identified uses of the substance or mixture and uses advised against

Laboratory reagent intended for in vitro diagnosis

1.3 Details of the supplier of the safety data sheet

Manufactured by:

DIMITRIADIS IOANNIS AND SON PC Production of Electrophoresis Products
Business Incubator Themi, Steliou Kazantzidi 47, 57001 THERMI PO BOX 8101,
GREECE Email: hellabio@hellabio.com. www.hellabio.com TEL (+30) 2311 999911

1.4 Emergency telephone number

Call (+30) 2311 999911 (Hellabio, Thessaloniki, Greece) for general information. In case of an emergency please contact your country's poison control center.

SECTION 2: Hazard identification

2.1 Classification of the substance or mixture according to Regulation (EC) 1272/2008

Human serum control

Signal word: Warning

| Health Hazards | Hazard Category |
|---|-----------------|
| Harmful if swallowed | Acute Tox. 4 |
| Contains biological material that may present a risk of infection | - |

2.2 Labeling Elements



SAFETY DATA SHEET

Date of issue: 15.01.2025

According to Regulation (EC) No. 2020/878 & 1272/2008

Page: 2/7

Human Antisera



Signal word: Danger

| Hazard statements | Hazard Class |
|-------------------|--------------|
| H302 | Acute Tox. 4 |

Precautions/Prevention

P203- Obtain, read, and follow all safety instructions before use

P280- Wear protective gloves, protective clothing, eye protection, face protection

P264- Wash hands thoroughly after handling

P270- Do not eat, drink, or smoke when using this product

2.3 Other Hazards

The serum has been tested for infectious agents; however it should still be treated as potentially infectious.

SECTION 3: Composition/information on ingredients

3.1 Substances / 3.2 Mixtures

Human control serum

| Component | Cas number | Conc.-% | Hazard Statement |
|-----------------|------------|---------|-------------------|
| Ethylene Glycol | 107-21-1 | 40% | H302- Acute Tox.4 |

All other components present no significant physical or chemical hazard at the concentration used.

SECTION 4. First aid measures

4.1 Description of first aid measures

Skin contact:

Remove contaminated clothing immediately. Rinse the affected skin or mucous membrane thoroughly for 15 minutes under running water.

Eye contact:

Rinse the affected eye(s) thoroughly under running water for at least 10 minutes. Further treatment by a physician should follow.

Inhalation:

Transfer the affected person(s) to fresh air. If symptoms persist, call a physician.



SAFETY DATA SHEET

Date of issue: 15.01.2025

According to Regulation (EC) No. 2020/878 & 1272/2008

Page: 3/7

Oral intake/ingestion:

Contact IMMEDIATELY with your country's poison control center

4.2 Most important symptoms and effects, both acute and delayed

In most cases, a rapid rinse with water is sufficient to alleviate the symptoms. However, prolonged exposure can lead to irritation and damage to the eyes and skin.

4.3 Indication of any immediate medical attention and special treatment needed

In case of oral intake/ingestion contact with country's poison control center and/or a physician IMMEDIATELY.

SECTION 5: Firefighting measures

5.1 Extinguishing media

All mixtures are non-flammable. Suitable extinguishing media: Water spray. Carbon dioxide (CO₂). Foam. Dry chemical. Not suitable extinguishing media: No data available.

5.2 Special hazards arising from the substance or mixture

No data available

5.3 Advice for firefighters

Wear protective suit with breathing apparatus

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Ensure adequate ventilation. Do not inhale vapors. Remove any contaminated clothing. Avoid contact with bare skin and use personal protective equipment (PPE). Clean with caution.

6.2 Environmental precautions

Avoid spilling into drains and waterways whenever possible.

6.3 Methods and material for containment and cleaning up

Wipe spillage with caution using absorbent paper. Do not touch with bare hands/skin. Rinse with water.

6.4 Reference to other sections

See section 3 for hazard information of each mixture

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Always follow the instructions for use. Each kit should be stored horizontally, as indicated on the label. Use personal protective equipment (Wear protective gloves, protective clothing, eye protection, face protection).



SAFETY DATA SHEET

Date of issue: 15.01.2025

According to Regulation (EC) No. 2020/878 & 1272/2008

Page: 4/7

7.2 Conditions for safe storage, including any incompatibilities

Store at -18-22°C for maximum 6 months.

7.3 Specific end use(s)

This kit is intended for in vitro diagnosis.

SECTION 8: Exposure controls/personal protection

8.1 Control parameters/ Occupational exposure limits (OELs)

| Component | Cas Number | ACGIH TLV® 8-hour TWA | CAL/OSHA PEL |
|-----------------|------------|--------------------------|----------------------------------|
| Ethylene Glycol | 107-21-1 | 10 mg/m ³ | 127 mg/m ³ (as vapor) |

8.2 Exposure controls

Use personal protective equipment (PPE)

Eye protection: Use safety glasses, avoid contact.Skin/Hand protection: Avoid direct contact of the product with skin, immediately remove any clothing soiled with the product and wash contaminated skin with soapy water. Use personal protective equipment, proper clothing and gloves.Respiratory protection: Avoid fumes, use in well-ventilated areas.Thermal hazards: Not applicable.

SECTION 9: Exposure controls/personal protection

9.1 Information on basic physical and chemical properties

Human serum control

| | |
|--|---------------------|
| Physical state | Liquid |
| Colour | Colourless |
| Odour | Odourless |
| Odour threshold | No data available |
| pH | - |
| Melting point/freezing point | No data available |
| Initial boiling point and boiling range | No data available |
| Flash point: | No data available |
| Evaporation rate: | No data available |
| Flammability (solid, gas) | Not flammable |
| Upper/lower flammability or explosive limits | No data available |
| Vapour pressure: | No data available |
| Vapour density: | No data available |
| Relative density: | No data available |
| Solubility(ies) | Miscible with water |
| Partition coefficient: n-octanol/water | No data available |
| Auto-ignition temperature | No data available |



SAFETY DATA SHEET

Date of issue: 15.01.2025

According to Regulation (EC) No. 2020/878 & 1272/2008

Page: 5/7

Decomposition temperature: No data available
Viscosity: No data available
Explosive properties: Not explosive

9.2 Other information

No other relevant information

SECTION 10: Stability and reactivity

10.1 Reactivity

The product is stable in conditions provided by the manufacturer.

10.2 Chemical stability

No decomposition if stored and applied as directed. The product is stable when normal handling in accordance with conditions provided by the manufacturer.

10.3 Possibility of hazardous reactions

No data available

10.4 Conditions to avoid

Avoid light and heat.

10.5 Incompatible materials

No data available

10.6 Hazardous decomposition products

No data available

SECTION 11: Toxicological information

There is no available data for the mixtures/buffers. Toxicity levels of chemicals used in mixtures are stated below. For Information on hazard classes as defined in Regulation (EC) No 1272/2008 see sections 2 and 3.

Acute toxicity

| Component | LD50 (Oral - Rat) |
|-----------------|-------------------|
| Ethylene glycol | 3000 mg/kg |

Skin corrosion/irritation: May irritate skin

Serious eye damage/irritation: Irritant for the eye

Respiratory or skin sensitization: No data available

Germ cell mutagenicity: No data available

Carcinogenicity: No data available

Reproductive toxicity: No data available

STOT-single exposure: No data available

STOT-repeated exposure: No data available

Aspiration hazard: No data available



SAFETY DATA SHEET

Date of issue: 15.01.2025

According to Regulation (EC) No. 2020/878 & 1272/2008

Page: 6/7

SECTION 12: Ecological information

12.1 Ecological Toxicity

Based on available data, the product does not present any significant risks to aquatic life or ecosystems under normal use and disposal conditions.

12.2 Persistence and degradability

No data available

12.3 Bioaccumulative potential

No data available

12.4 Mobility in soil

No data available

12.5 Results of PBT and vPvB assessment

No data available

12.6 Endocrine disrupting properties

No data available

12.7 Other adverse effects

No other relevant information

SECTION 13: Disposal considerations

13.1 Waste treatment methods

It is crucial to handle and dispose of chemicals in compliance with local, regional, and national regulations to prevent environmental contamination or harm to public health. Users are advised to exercise caution when managing waste, ensuring that disposal methods minimize risks to human health and the environment. Proper waste treatment can involve neutralization, incineration, or disposal at a licensed facility designed to handle hazardous substances. Always consult local waste disposal guidelines and regulatory authorities to ensure full compliance with your country's laws regarding hazardous material disposal.

SECTION 14: Transport information

14.1 UN number or ID number

Not applicable

14.2 UN proper shipping name

Not applicable

14.3 Transport hazard class(es)

Not applicable

14.4 Packing group

No limits

14.5 Environmental hazards

Avoid spillage to the environment

14.6 Special precautions for user

Not applicable

14.7 Maritime transport in bulk according to IMO instruments

Not applicable



SAFETY DATA SHEET

Date of issue: 15.01.2025

According to Regulation (EC) No. 2020/878 & 1272/2008

Page: 7/7

SECTION 15: Regulatory information

This Material Safety Data Sheet (MSDS) was prepared in accordance with Regulation (EC) No 1907/2006 (REACH), No 1272/2008 and No 2020/878.

The supplier has not carried out a chemical evaluation or safety assessment for the substance or mixture.

SECTION 16: Other information

The information provided in this Material Safety Data Sheet (MSDS) is intended solely for safety guidance and is based on the best available knowledge at the time of publication. It is not intended to be a warranty. The manufacturer or distributor does not assume any liability for the misuse or improper handling of the product. Always refer to the latest regulatory information and consult with safety experts when necessary.