



## HELLABIO PRECAST GELS

### Separation Ability (MW size range in bp) of Gels

The Separation ability depends on the quality and the concentration of the Agarose

Cod	Size range in bp
D01	10-300
D02	20-500
D03	50-100
D04	50-2000
D05	100-3500
D06	200-6000
D07	300-8000
D08	500-12 000
D09	1000-20000

**Tabel 2**  
**Thickness of the Gels: 4 mm. Free of Ethidium Bromide\*\***

Subcode	Width/Length of the Gels (mm)	Lot	Distance between wells	Number of Teeth's	Number of Samples	Width/Length of the wells	Price in €**
G01	59/83	1	2	9	9	1x4	55
G01	59/83	2	3	8	8	1x4	54
G01	59/83	3	4	8	8	1x4	54
G01	59/83	4	2	11	11	1x3	57
G01	59/83	5	3	9	9	1x3	55
G02	83/59	1	2	13	13	1x4	59
G02	83/59	2	3	11	11	1x4	57
G02	83/59	3	4	10	10	1x4	56
G02	83/59	4	2	15	15	1x3	61
G02	83/59	5	3	13	13	1x3	59
G03	70/90	1	2	11	11	1x4	61
G03	70/90	2	3	9	9	1x4	58
G03	70/90	3	4	8	8	1x4	57
G03	70/90	4	2	13	13	1x3	62
G03	70/90	5	3	11	11	1x3	60
G04	90/70	1	2	14	14	1x4	63
G04	90/70	2	3	12	12	1x4	61
G04	90/70	3	4	10	10	1x4	59
G04	90/70	4	2	17	17	1x4	68
G04	90/70	5	3	14	14	1x4	63

\*\* Plus 5% with Ethidium

*How to order ( example):*

<i>Cod</i>	<i>Subcode</i>	<i>Lot</i>	<i>Ethidium Bromide</i>	<i>Price /kit</i>
<i>D01</i>	<i>G01</i>	<i>5</i>	<i>no</i>	<i>55</i>
<i>D01</i>	<i>G01</i>	<i>5</i>	<i>yes</i>	<i>57.75</i>