



**PRICE LIST
OF HELLABIO ELECTROPHORESIS KITS AND EQUIPMENTS**

IMMUNOFIXATION KITS

COD:	LOT	TESTS/ KIT	CONC. BUFFER FOR LITER	INTENDED TO BE USED
IFE 01	01	10	1	Immunofixation Electrophoresis of Proteins
	02	10	3	
IFE D01	03	20	1	
	04	20	3	

PROTEIN ELECTROPHORESIS KITS

COD:	LOT	TESTS/ KIT	CONC. BUFFER FOR LITER	COD: OF TEMPLATES	INTENDED TO BE USED
PE10	21	100	1	E10	Electrophoresis of Proteins
	22	100	3	E10	
	23	100	1	PE011	
	24	100	3	PE011	
	25	100	1	T01	
	26	100	3	T01	
PE14	27	140	1	01.15	
PE14	28	140	3	01.15	

HEMOGLOBIN ELECTROPHORESIS KITS

COD:	LOT:	TESTS/ KIT	CONC. BUFFER FOR LITER	COD: OF TEMPLATES	INTENDED TO BE USED
HE10	51	100	1	E10	Electrophoresis of
	52	100	3	E10	
	53	100	1	PE011	
	54	100	3	PE011	
	55	100	1	T01	
	56	100	3	T01	
HE 14	57	140	1	01.15	
	58	140	3	01.15	
HE10	71	100	1	E10	Electrophoresis of
	72	100	3	E10	
	73	100	1	PE011	
	74	100	3	PE011	
	75	100	1	T01	
	76	100	3	T01	
HE 14	57	140	1	01.15	
	58	140	3	01.15	

LIPOPROTEIN ELECTROPHORESIS KITS

COD:	LOT:	TESTS/KIT	CONC. BUFFER FOR LITER	COD: OF TEMPLATES	INTENDED TO BE USED
LE10	91	100	1	E10	Lipoprotein electrophoresis
	92	100	3	E10	
	93	100	1	PE011	
	94	100	3	PE011	
	95	100	1	T01	
	96	100	3	T01	
LE14	97	140	1	01.15	
	98	140	3	01.15	

IMMUNOELECTROPHORESIS KITS

COD:	LOT:	TESTS/KIT	CONC. BUFFER FOR LITER	COD OF TEMPLATES	INTENDED TO BE USED
IE 01	111	10	1	IE1	Immunoelectrophoresis of proteins
	112	10	3	IE1	

SERUM AND HEMOGLOBIN CONTROLS

Cod	Lot	Name of the Product	Package
HS01	150	SERUM CONTROL 0.5ML	1
HP01	151	PLASMA CONTROL 0.5ML	1
HbS01	152	HEMOGLOBIN S (HbS) 0.25 ML	1
HbF 01	153	HEMOGLOBIN F(HbF) 0.25 ML	1
•HbA201	154	HEMOGLOBIN A2 (HbA2) 0.25 ML	1
HbA 01	155	HEMOGLOBIN A (HbA) NORMAL 0.25ML	1

CRYOGLOBULIN KIT

COD:	LOT:	TESTS/KIT	INTENDED TO BE USED
HCG 01	01	108	Detection of Cryoglobulins

KIT COMPONENTS IN SEPARATED FORM

Cod	Accessories	Exp. Date
EA01	Conc. Protein Electrophoresis Buffer for 1 liter	24 months
EA02	Conc. Protein Electrophoresis Buffer for 3 liter	24 months
EA03	Conc. Hemoglobin Electrophoresis Buffer for 1 liter	24 months
EA04	Conc. Hemoglobin Electrophoresis Buffer for 3 liter	24 months
EA05	Gel Blotter Strips (20 pieces)	-
EA06	Templates (10 pieces)	-
EA07	Conc. Staining Solution for 300 ml	24 months
EA08	Protein diluent solution 5 ml	24 months
EA09	10 Agarose Gels Cod PE 10	15 months
EA10	10 Agarose Gels Cod PE 14	15 months
EA11	10 Agarose Gels Cod HE 10	15 months
EA12	10 Agarose Gels Cod HE 14	15 months
EA13	Haemolysing Solution 35 ml	24 months

ANIMAL BLOODS

Cod	Name of the Products	Package
SBS 01	Sheep Blood in Alsever	50 ml
SBS 01	Sheep Blood sterile	50 ml
HB 01	Horse Blood sterile	350 ml
HBH 01	Horse Blood Hemolyzed sterile	350 ml

ANTISERA

Cod	Name of the Product	Properties	Package
SAHG2	Sheep anti Human IgG (H+L chain)	IgG grade	2ml
SAHG10	Sheep anti Human IgG (H+L chain)	IgG grade	10ml
SAHG50	Sheep anti Human IgG (H+L chain)	IgG grade	50ml
SAHW2	Sheep anti Human Whole serum	IgG grade	2 ml
SAHW10	Sheep anti Human Whole serum	IgG grade	10 ml
SAHW50	Sheep anti Human Whole serum	IgG grade	50 ml
SARG2	Sheep anti Rabbit IgG (H+L)	IgG grade	2 ml
SARG10	Sheep anti Rabbit IgG (H+L)	IgG grade	10 ml
SARG50	Sheep anti Rabbit IgG (H+L)	IgG grade	50 ml

PLASTIC TRAYS

Cod	Name of the Product	Package
PTCG 01	Clear plastic tray for Ouchterlony	100 plates

HELLABIO ELECTROPHORESIS EQUIPMENTS

A. HORIZONTAL ELECTROPHORESIS SYSTEM

1. HORIZONTAL ELECTROPHORESIS TANK



Model	Cod No	Lot	Name of the Product/Description	Package
HET01	ET01	01	Electrophoresis tank, one bridge , wire	1
HET01	ET01	02	Electrophoresis tank, one bridge	1
HET02	DET01	01	Electrophoresis tank with wire, two bridges	1


2. Accessories for Horizontal Electrophoresis Tank

Cod	Lot	Name of the Product/Description	Package
SB01	01	Five Staining-Distaining Baths +Basis	1
FT01	01	Film Transporter	2

Cod	Lot	Name of the Product	Package
HSCANSOI	01	HellabioScan Software (without Hardware)	1
HPS01	01	Power Supply 25-200 volts	1
HFD	01	Del Drier	1
HPS500	01	Power Supply 30-500 volts	1
HTC	01	TANK COOLER	1

		
Power Supply	Gel Drier	TANK COOLER. Store it at -8°C. Before start running electrophoresis put the electrophoresis tank on the cooler box.

**4. Dams (color: Red /Black.
Thickness: 4mm)**



Model	Cod	Lot	Name of the product	Size Height/ width (mm)	Package
HD	01	01	Dam	40X59	2
HD	01	02	Dam	40X83	2
HD	01	03	Dam	40X70	2
HD	01	04	Dam	40X90	2

5. UV Transparent Gel Carrier

Model	Cod	Lot	Name of the product	Size Width/length (mm)	Package
HUVGCOI	GC01	01	UV Transparent Gel Carrier	59X83	1
HUVGCOI	GC01	02	UV Transparent Gel Carrier	83X59	1
HUVGCOI	GC01	03	UV Transparent Gel Carrier	70X90	1
HUVGCOI	GC01	04	UV Transparent Gel Carrier	90X70	1

6. Wedge



Model	Cod	Lot	Name of the product	Size Height/ length (mm)	Package
HW	W01	01	Wedge	24X59	1
HW	W01	02	Wedge	24X83	1
HW	W01	03	Wedge	24X70	1
HW	W01	04	Wedge	24X90	1

C. VERTICAL ELECTROPHORESIS SYSTEM

Model	Cat No	Lot	Name of the Product	Gel size	Package
HVES01	01	01	Vertical Electrophoresis Tank, Stand to Precast gels, Wells Forming Comb Strips for precasting the Gel. Glass plates to precast gels	15X17.3	1
HVES02	01	02	Vertical Electrophoresis Tank, Stand to Precast gels, Wells Forming Comb, Strips for precasting the Gel. Glass plates to precast gels	11X13.5	1
HVES01	02	01	Stand to precast gel		1
HVES01	03	01	Glass plates		2
HVES01	04	01	Strips to precast gels		3



Model	Cod	Lot	Name of the product	Package
HDRB	DRB01	01	Dark room Box	1

HELLABIO PRECAST GELS

Separation Ability (MW size range in bp) of Gels

The Separation ability depends on the quality and the concentration of the Agarose

Cod	Size range in bp
D01	10-300
D02	20-500
D03	50-100
D04	50-2000
D05	100-3500
D06	200-6000
D07	300-8000
D08	500-12 000
D09	1000-20000

Hellabio Precast Gels for submarine Electrophoresis [Thickness of the Gels: 4 mm. Free of Ethidium Bromide

Sub code	Lot	Width/Length of the Gels (mm)	Width/Length of the wells (mm)	Distance between wells (mm)	Number of teeth
G01	1	59/83	1x4	2	9
	2	59/83	1x4	3	8
	3	59/83	1x4	4	8
	4	59/83	1x3	2	11
	5	59/83	1x3	3	9
G02	1	83/59	1x4	2	13
	2	83/59	1x4	3	11
	3	83/59	1x4	4	10
	4	83/59	1x3	2	15
	5	83/59	1x3	3	13
G03	1	70/90	1x4	2	11
	2	70/90	1x4	3	9
	3	70/90	1x4	4	8
	4	70/90	1x3	2	13
	5	70/90	1x3	3	11
G04	1	90/70	1x4	2	14
	2	90/70	1x4	3	12
	3	90/70	1x4	4	10
	4	90/70	1x4	2	17
	5	90/70	1x4	3	14

How to order Precast Gels (Example):

Code	Sub code	Lot	Ethidium Bromide
D01	G01	5	no
D01	G01	5	yes

--	--