



# SAFETY DATA SHEET

Date of issue: 15.01.2025

According to Regulation (EC) No. 2020/878 & 1272/2008

Page: 1/7

## SECTION 1: Identification of the substance/mixture and of the company

### 1.1 Product Identifier

REF: HSS-FR

Product name: FAT RED STAINING SOLUTION

REACH registration number: The annual tonnage does not require registration

### 1.2 Relevant identified uses of the substance or mixture and uses advised against

Laboratory reagents intended to assist with in vitro diagnosis

### 1.3 Details of the supplier of the safety data sheet

Manufactured by:

DIMITRIADIS IOANNIS AND SON PC Production of Electrophoresis Products  
Business Incubator Thermi, Steliou Kazantzidi 47, 57001 THERMI PO BOX 8101,  
GREECE Email: hellabio@hellabio.com. www.hellabio.com TEL (+30) 2311 999911

### 1.4 Emergency telephone number

Call (+30) 2311 999911 (Hellabio, Thessaloniki, Greece) for general information. In case of an emergency please contact your country's poison control center.

## SECTION 2: Hazard identification

### 2.1 Classification of the substance or mixture according to Regulation (EC) 1272/2008

Staining solution

**Signal word: Danger**

Health Hazards	Hazard Category
Harmful if swallowed	Acute Tox. 4
Highly Flammable liquid and vapor	Flam. Liq. 2
Toxic if swallowed	Acute Tox. 3
Toxic in contact with skin	Acute Tox. 3
Toxic if inhaled	Acute Tox. 3
Causes damage to organs	STOT SE 1

### 2.2 Labeling Elements

# SAFETY DATA SHEET

Date of issue: 15.01.2025

According to Regulation (EC) No. 2020/878 & 1272/2008

Page: 2/7

## Staining solution



Signal word: Danger

Hazard statements	Hazard Class
H225	Flam. Liq. 2
H301	Acute Tox. 3
H303	Harmful if swallowed
H311	Acute Tox. 3
H331	Acute Tox. 3
H370	STOT SE 1

Precautions/Prevention
P203- Obtain, read, and follow all safety instructions before use
P210- Keep away from heat, hot surface, sparks, open flames and other ignition sources. No smoking.
P233- Keep container tightly closed.
P280- Wear protective gloves, protective clothing, eye protection, face protection
P264- Wash hands thoroughly after handling
P261- Avoid breathing dust/fume/gas/mist/vapours/spray
P270- Do not eat, drink or smoke when using this product.
P304+P340: If inhaled, remove person to fresh air and keep comfortable for breathing.

## 2.3 Other Hazards

Not applicable

## SECTION 3: Composition/information on ingredients

### 3.1 Substances / 3.2 Mixtures

#### Staining solution

Component	Cas number	Conc.-%	Hazard Statement
Methanol	77-92-9	95-99%	H225- Highly Flammable liquid and vapor H301- Toxic if swallowed H311- Toxic in contact with skin H331- Toxic if inhaled H370-Causes damage to organs



# SAFETY DATA SHEET

Date of issue: 15.01.2025

According to Regulation (EC) No. 2020/878 & 1272/2008

Page: 3/7

All other components present no significant physical or chemical hazard at the concentration used.

## SECTION 4: First aid measures

### 4.1 Description of first aid measures

#### Skin contact:

Remove contaminated clothing immediately. Rinse the affected skin or mucous membrane thoroughly for 15 minutes under running water.

#### Eye contact:

Rinse the affected eye(s) thoroughly under running water for at least 10 minutes. Further treatment by a physician should follow.

#### Inhalation:

Transfer the affected person(s) to fresh air. If symptoms persist, call a physician.

#### Oral intake/ingestion:

Contact IMMEDIATELY with your country's poison control center

### 4.2 Most important symptoms and effects, both acute and delayed

In most cases, a rapid rinse with water is sufficient to alleviate the symptoms. However, prolonged exposure can lead to irritation and damage to the eyes and skin.

### 4.3 Indication of any immediate medical attention and special treatment needed

In case of oral intake/ingestion contact with country's poison control center and/or a physician IMMEDIATELY.

## SECTION 5: Firefighting measures

### 5.1 Extinguishing media

Staining solution is highly flammable. Suitable extinguishing media: Carbon dioxide (CO<sub>2</sub>). Foam. Dry chemical. Not suitable extinguishing media: Water.

### 5.2 Special hazards arising from the substance or mixture

No data available

### 5.3 Advice for firefighters

Wear protective suit with breathing apparatus

## SECTION 6: Accidental release measures

### 6.1 Personal precautions, protective equipment and emergency procedures

Ensure adequate ventilation. Do not inhale vapors. Remove any contaminated clothing. Avoid contact with bare skin and use personal protective equipment (PPE). Clean with caution.

### 6.2 Environmental precautions



# SAFETY DATA SHEET

Date of issue: 15.01.2025

According to Regulation (EC) No. 2020/878 & 1272/2008

Page: 4/7

Avoid spilling into drains and waterways whenever possible. Buffers containing sodium azide pose a danger to the environment.

## 6.3 Methods and material for containment and cleaning up

Wipe spillage with caution using absorbent paper. Do not touch with bare hands/skin. Rinse with water.

## 6.4 Reference to other sections

See section 3 for hazard information of each mixture

## SECTION 7: Handling and storage

### 7.1 Precautions for safe handling

Always follow the instructions for use. Each kit should be stored horizontally, as indicated on the label. Use personal protective equipment (Wear protective gloves, protective clothing, eye protection, face protection).

### 7.2 Conditions for safe storage, including any incompatibilities

Store in a dry, cool and well-ventilated area between 15-25°C.

### 7.3 Specific end use(s)

This kit is intended for in vitro diagnosis

## SECTION 8: Exposure controls/personal protection

### 8.1 Control parameters/ Occupational exposure limits (OELs)

Component	Cas Number	ACGIH TLV <sup>©</sup> 8-hour TWA	CAL/OSHA PEL
Methanol	67-56-1	262 mg/m <sup>3</sup>	260 mg/m <sup>3</sup>

### 8.2 Exposure controls

Use personal protective equipment (PPE)

Eye protection: Use safety glasses, avoid contact.

Skin/Hand protection: Avoid direct contact of the product with skin, immediately remove any clothing soiled with the product and wash contaminated skin with soapy water. Use personal protective equipment, proper clothing and gloves.

Respiratory protection: Avoid fumes, use in well-ventilated areas.

Thermal hazards: Keep away from fire and heat surfaces.

## SECTION 9: Exposure controls/personal protection

### 9.1 Information on basic physical and chemical properties

#### Staining solution

Physical state

Liquid

Colour

Red



# SAFETY DATA SHEET

Date of issue: 15.01.2025

According to Regulation (EC) No. 2020/878 & 1272/2008

Page: 5/7

Odour	light
Odour threshold	No data available
pH	<2
Melting point/freezing point	No data available
Initial boiling point and boiling range	No data available
Flash point:	No data available
Evaporation rate:	No data available
Flammability (solid, gas)	Highly Flammable
Upper/lower flammability or explosive limits	No data available
Vapour pressure:	No data available
Vapour density:	No data available
Relative density:	No data available
Solubility(ies)	Miscible with water
Partition coefficient: n-octanol/water	No data available
Auto-ignition temperature	No data available
Decomposition temperature:	No data available
Viscosity:	No data available
Explosive properties:	Not explosive

## 9.2 Other information

No other relevant information

## SECTION 10: Stability and reactivity

### 10.1 Reactivity

The product is stable in the conditions provided by the manufacturer.

### 10.2 Chemical stability

No decomposition if stored and applied as directed. The product is stable when normal handling is in accordance with conditions provided by the manufacturer.

### 10.3 Possibility of hazardous reactions

No data available

### 10.4 Conditions to avoid

Avoid light and heat.

### 10.5 Incompatible materials

No data available

### 10.6 Hazardous decomposition products

No data available

## SECTION 11: Toxicological information



# SAFETY DATA SHEET

Date of issue: 15.01.2025

According to Regulation (EC) No. 2020/878 &amp; 1272/2008

Page: 6/7

There is no available data for the mixtures/buffers. Toxicity levels of chemicals used in mixtures are stated below. For Information on hazard classes as defined in Regulation (EC) No 1272/2008 see sections 2 and 3.

## Acute toxicity

Component	LD50 (Oral - Rat)
Methanol	>5000 mg/kg

Skin corrosion/irritation: May irritate skin  
Serious eye damage/irritation: Irritant for the eye  
Respiratory or skin sensitization: No data available  
Germ cell mutagenicity: No data available  
Carcinogenicity: No data available  
Reproductive toxicity: No data available  
STOT-single exposure: No data available  
STOT-repeated exposure: No data available  
Aspiration hazard: No data available

## SECTION 12: Ecological information

### 12.1 Ecological Toxicity

No data available.

### 12.2 Persistence and degradability

No data available

### 12.3 Bioaccumulative potential

No data available

### 12.4 Mobility in soil

No data available

### 12.5 Results of PBT and vPvB assessment

No data available

### 12.6 Endocrine disrupting properties

No data available

### 12.7 Other adverse effects

No other relevant information

## SECTION 13: Disposal considerations

### 13.1 Waste treatment methods

It is crucial to handle and dispose of chemicals in compliance with local, regional, and national regulations to prevent environmental contamination or harm to public health. Users are advised to exercise caution when managing waste, ensuring that disposal methods minimize risks to human health and the environment. Proper waste treatment can involve neutralization, incineration, or disposal at a licensed facility designed to handle hazardous substances. Always consult local waste disposal guidelines and regulatory authorities to ensure full compliance with your country's laws regarding hazardous material disposal.



# SAFETY DATA SHEET

Date of issue: 15.01.2025

According to Regulation (EC) No. 2020/878 & 1272/2008

Page: 7/7

## SECTION 14: Transport information

### 14.1 UN number or ID number

Not applicable

### 14.2 UN proper shipping name

Not applicable

### 14.3 Transport hazard class(es)

Not applicable

### 14.4 Packing group

No limits

### 14.5 Environmental hazards

Avoid spillage to the environment

### 14.6 Special precautions for user

Not applicable

### 14.7 Maritime transport in bulk according to IMO instruments

Not applicable

## SECTION 15: Regulatory information

This Material Safety Data Sheet (MSDS) was prepared in accordance with Regulation (EC) No 1907/2006 (REACH), No 1272/2008 and No 2020/878.

The supplier has not carried out a chemical evaluation a safety assessment for the substance or mixture.

## SECTION 16: Other information

The information provided in this Material Safety Data Sheet (MSDS) is intended solely for safety guidance and is based on the best available knowledge at the time of publication. It is not intended to be a warranty. The manufacturer or distributor does not assume any liability for the misuse or improper handling of the product. Always refer to the latest regulatory information and consult with safety experts when necessary.